


Product Name: Recombinant Human EFNA3 Protein
Catalog #: PEH5722



Summary

Name	EFNA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P52797
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGPGGALGNR HAVYWNSSNQ HLRREGYTVQ VNVNDYLDIY CPHYNSSGVG PGAGPGPGGG AEQYVLYMVS RNGYRTCNAS QGFKRWECNR PHAPHSPIKF SEKFQRYSAF SLGYEFHAGH EYYYISTPTH NLHWKCLRMK VFVCCASTSH SGEKPVP TLP QFTMGPNVKI NVLEDFEGEN PQVPKLEKSI SG

image



Background

Alternative Names	Ephrin-A3 (EFL-2) (EHK1 ligand) (EHK1-L) (EPH-related receptor tyrosine kinase ligand 3) (LERK-3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EFNA3 Protein
Catalog #: PEH5722

