


Product Name: Recombinant Human CADM2 Protein
Catalog #: PEH5730



Summary

Name	CADM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3J6
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFSSPVMEGD LMQLTCKTSG SKPAADIRWF KNDKEIKDVK YLKEEDANRK TFTVSSTLDF RVDRSDDGVA VICRVDHESL

image



Background

Alternative Names	Cell adhesion molecule 2 (Immunoglobulin superfamily member 4D) (IgSF4D) (Nectin-like protein 3) (NECL-3) (Synaptic cell adhesion molecule 2) (SynCAM 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use