


Product Name: Recombinant Human DDX56 Protein
Catalog #: PEH5736



Summary

Name	DDX56_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NY93
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QESQLPGPDQ LQQFQVVCET EEDKFLLLYA LLKLSLIRGK SLLFVNTLER SYRLRLFLEQ FSIPTCVLNG ELPLRSRCHI ISQFNQGFYD CVIATDAEVL GAPVKGKRRG RGPKGDKASD PEAGVARGID FHHVSAVLNF DLPPTPEAYI HRAGRTARAN NPGIVLTFVL PTEQFHLGKI EELLSGENRG PILLPYQFRM EEIEGFRYRC RDAMRSVTKQ AIREARLKEI KEELLHSEKL KTYFEDNPRD

image



Background

Alternative Names	Probable ATP-dependent RNA helicase DDX56 (EC 3.6.4.13) (ATP-dependent 61 kDa nucleolar RNA helicase) (DEAD box protein 21) (DEAD box protein 56)
--------------------------	---

Note

Product Name: Recombinant Human DDX56 Protein
Catalog #: PEH5736



For Research Use Only , Not for Diagnostic Use