


Product Name: Recombinant Human KLK3 Protein
Catalog #: PEH5741



Summary

Name	KLK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	P07288
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CEKHSQPWQV LVASRGRAVC GGVLVHPQWV LTAAH CIRNK SVILLGRHSL FHPEDTGQVF QVSHSFPHPL YDMSLLKNRF LRPGDDSSHD LMLLRLESPA ELTDAVKVMD LPTQEPALGT TCYASGWGSI EPEEFLTPKK LQCVDLHVIS NDVCAQVHPQ KVTKFMCLCAG RWTGGKSTCS GDSGGPLVCN GVLQGITSWG

image



Background

Alternative Names Prostate-specific antigen (PSA) (EC 3.4.21.77) (Gamma-seminoprotein) (Seminin) (Kallikrein-3) (P-30 antigen) (Semenogelase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KLK3 Protein
Catalog #: PEH5741

