


Product Name: Recombinant Human GPAT2 Protein
Catalog #: PEH5753



Summary

Name	GPAT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NUI2
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRSLQHSLS LLRAHVALLR IRQGDLLVVP QPGPGLTHLA QLSAELLPVF LSEAVGACAV RGLLAGRVPP QGPWELQGIL LLSQNELYRQ ILLLMHLLPQ DLLLLKPCQS SYCYCQEVLD RLIQCGLLVA EETPGSRPAC DTGRQRLSRK LLWKPSGDFT DSDSDDFGEA DGRYFRLSQQ SHCPDFFLFL CRLLSPLLKA FAQAAAFRLRQ GQLPDTELG YTEQLFQFLQA TAQEEGIFEC

image



Background

Alternative Names	Glycerol-3-phosphate acyltransferase 2, mitochondrial (GPAT-2) (EC 2.3.1.15) (1-acylglycerol-3-phosphate O-acyltransferase GPAT2) (EC 2.3.1.51) (xGPAT1)
--------------------------	--

Product Name: Recombinant Human GPAT2 Protein
Catalog #: PEH5753



Note

For Research Use Only , Not for Diagnostic Use