


Product Name: Recombinant Human GZMA Protein
Catalog #: PEH5757




Summary

Name	GRAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P12544
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IIGGNEVTPH SRPYMVLLSL DRKTICAGAL IAKDWVLTAH HCNLNKRSQV ILGAHSITRE EPTKQIMLVK KEFPYPCYDP ATREGDLKLL QLMEKAKINK YVTILHLPKK GDDVKPGTMC QVAGWGRTHN SASWSDTLRE VNITIIDRKV CNDRNHYNFN PVIGMNMVCA GSLRGGRDSC NGDSGSPLLC EGVFRGVTSF GLENKCGDPR GPGVYILLSK KHLNWIIMTI KGAV

image



SDS



WB

Background

Alternative Names	Granzyme A (EC 3.4.21.78) (CTL tryptase) (Cytotoxic T-lymphocyte proteinase 1) (Fragmentin-1) (Granzyme-1) (Hanukkah factor) (H factor) (HF)
--------------------------	--

Note

Product Name: Recombinant Human GZMA Protein
Catalog #: PEH5757



For Research Use Only , Not for Diagnostic Use