


Product Name: Recombinant Human HROB Protein
Catalog #: PEH5758



Summary

Name	HROB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3J3
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRSLEDIMVS APQTPTHGAL AKFQTEIVAS SQASVEEDFG RGPWLTMKST LGLDERDPSC FLCTYSIVMV LRKQAALKQL PRNKVPNMAV MIKSLTRSTM DASVVFKDPT GEMQGTVHRL LLETQNELK PGSVLLKQI GVFPSPSLRNH YLNVTNNLV HIYSPDSGDG SFLKPSQFPF KDSGSFQHDV AAKPEEGFRT AQNLEAEASP EEELPEADDL DGLLSELPED

image



Background

Alternative Names Homologous recombination OB-fold protein

Note

Product Name: Recombinant Human HROB Protein
Catalog #: PEH5758



For Research Use Only , Not for Diagnostic Use