


Product Name: Recombinant Human CSTF2T Protein
Catalog #: PEH5770



Summary

Name	CSTFT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0L4
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDRSLRSVFV GNIPYEATEE QLKDIFSEVG SVVSFRLVYD RETGKPKGYG FCEYQDQETA LSAMRNLNGR EFSGRALRVD NAASEKNKEE LKSLGPAAPI IDSPYGDPID PEDAPESITR AVASLPPEQM FELMKQMKLC VQNSHQEARN MLLQNPQLAY ALLQAQVVMR IMDPEIALKI LHRKIHVTPV IPGKSQSVSV SGPGPGPGPG LCPGPNVLLN QQNPPAPQPQ HLARRPVKDI PPLMQTPIQG GIPAPGPIPA

image



Background

Alternative Names	Cleavage stimulation factor subunit 2 tau variant (CF-1 64 kDa subunit tau variant) (Cleavage stimulation factor 64 kDa subunit tau variant) (CSTF 64 kDa subunit tau variant) (TauCstF-64)
--------------------------	---

Product Name: Recombinant Human CSTF2T Protein
Catalog #: PEH5770



Note

For Research Use Only , Not for Diagnostic Use