



Product Name: Recombinant Human MCOLN2 Protein
Catalog #: PEH5777



Summary

Name	MCLN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZK6
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IFFAINQYHQ LKDLTLGLG YGENEDNRIG LKVCQHYKK GTMFPSNETL NIDNDVELDC VQLDLQDLSK KPPDWKNSSF FRLEFYRLQ VEISFHLKGI DLQTIHSREL PDCYVFQNTI IFDNKAHSGK IKIYFSDAK IEECKDLNIF

image



Background

Alternative Names Mucolipin-2 (Transient receptor potential channel mucolipin 2) (TRPML2)

Note

For Research Use Only , Not for Diagnostic Use