


Product Name: Recombinant Human CPHXL Protein
Catalog #: PEH5795



Summary

Name	CPHXL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A0A1W2PPM1
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAEEDHHNEE RQTKNKRKTK HRHKFSEELL QELKEIFGEN CYPDYTTTRKT LAIKFDCPVN VIDNWFQNK RARLPPAERRR IFVLQKKHDF PVQAHSFLSC QETQAAAHNY ATKQSLSGAQ RALMRRAGCS HLEKQWIPSQ EMGYNCFSLE NQETPSQQVG PQCSYLEKPG IPSQQVGSQC SYLEKLIGPS QQVASQSSYL VTGTEKHPGC AMGYGGDTGS GHSGSGHSTA YHFLSYNSAE CLHPPPSSVP YFHGERTETK

image



Background

Alternative Names	Cytoplasmic polyadenylated homeobox-like (Cytoplasmic polyadenylated homeobox 1)
--------------------------	--

Product Name: Recombinant Human CPHXL Protein
Catalog #: PEH5795



Note

For Research Use Only , Not for Diagnostic Use