


Product Name: Recombinant Human EFNA2 Protein
Catalog #: PEH5796



Summary

Name	EFNA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43921
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAEDAARANS DRYAVYWNRS NPRFHAGAGD DGGGYTVEVS INDYLDIYCP HYGAPLPPAE RMEHYVLYMV NGEGHASCDH RQRGFKRWEC NRPAAPGGPL KFSEKFQLFT PFSLGFEFRP GHEYYYISAT PPNVAVDRPCL RLKVYVRPTN ETLYEAPEPI FTSN

image



Background

Alternative Names Ephrin-A2 (EPH-related receptor tyrosine kinase ligand 6) (LERK-6) (HEK7 ligand) (HEK7-L)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EFNA2 Protein
Catalog #: PEH5796

