


Product Name: Recombinant Human LYPD1 Protein
Catalog #: PEH5804



Summary

Name	LYPD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8N2G4
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQIQCYQCEE FQLNNDCCSP EFIVNCTVNV QDMCQKEVME QSAGIMYRKS CASSAACLIA SAGYQSFCSG GKLNSVCISC CNTPLCNGPR PKKRGS

image



Background

Alternative Names	Ly6/PLAUR domain-containing protein 1 (Putative HeLa tumor suppressor) (PHTS)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use