


Product Name: Recombinant Human CATSPER3 Protein
Catalog #: PEH5808



Summary

Name	CTSR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XQ3
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEDSIRKFER ELMLEQQEML MGEKQVILQR QQEEISRLMH IQKNADCTSF SELVENFKKT LSHTDPMVLD DFGTSLPFID IYFSTLDYQD TTVHKLQELY YEIVHVLSLM LEDLPQEKPQ

image



Background

Alternative Names Cation channel sperm-associated protein 3 (CatSper3) (Ca(v)-like protein)
(One-repeat calcium channel-like protein)

Note

For Research Use Only , Not for Diagnostic Use