



Product Name: Recombinant Human CTSW Protein
Catalog #: PEH5811



Summary

Name	CATW_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P56202
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPFSCDWRKV ASAISPIKDQ KNCNCCWAMA AAGNIETLWR ISFWDFVDVS VQELDCGRC GDGCHGGFVW DAFITVLNNS GLASEKDYPF QGKVRHRCH PKKYQKVAWI QDFIMLQNE HRIAQYLATY GPITVTINMK PLQLYRKGVI KATPTTCDPQ LVDHSVLLVG FGSVKSEEGI WAETVSSQSQ PQPPHPTPYW ILKNSWGAQW GEKGYFRLHR GSNTCGITKF PLTARVQKPD MKPRVSCPP

image



Background

Alternative Names Cathepsin W (EC 3.4.22.-) (Lymphopain)

Note

Product Name: Recombinant Human CTSW Protein
Catalog #: PEH5811



For Research Use Only , Not for Diagnostic Use