


Product Name: Recombinant Human CCL26 Protein
Catalog #: PEH5814



Summary

Name	CCL26_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y258
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRGSDISKTC CFQYSHKPLP WTWVRSYEFT SNSCSQRAVI FTTKRGKKVC THPRKKWVQK YISLLKTPKQ L

image



Background

Alternative Names	C-C motif chemokine 26 (CC chemokine IMAC) (Eotaxin-3) (Macrophage inflammatory protein 4-alpha) (MIP-4-alpha) (Small-inducible cytokine A26) (Thymic stroma chemokine-1) (TSC-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use