


Product Name: Recombinant Human CCL27 Protein
Catalog #: PEH5818



Summary

Name	CCL27_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9Y4X3
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLLPPSTACC TQLYRKPLSD KLLRKVIQVE LQADGDCHL QAFVLHLAQR SICIHPQNPS LSQWFEHQER KLHGTLPKLN FGMLRKMG

image



Background

Alternative Names	C-C motif chemokine 27 (CC chemokine ILC) (Cutaneous T-cell-attracting chemokine) (CTACK) (ESkine) (IL-11 R-alpha-locus chemokine) (Skinkine) (Small-inducible cytokine A27)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use