


Product Name: Recombinant Human NETO2 Protein
Catalog #: PEH5826



Summary

Name	NETO2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NC67
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGQNIGIKHI PATQCGIWVR TSNGGHFASP NYPDSYPPNK ECIYILEAAP RQRIELTFDE HYYIEPSFEC RFDHLEVRDG PFGFSPLIDR YCGVKSPPLI RSTGRFMWIK FSSDEELEG L GFRAKYSFIP

image



Background

Alternative Names	Neuropilin and tolloid-like protein 2 (Brain-specific transmembrane protein containing 2 CUB and 1 LDL-receptor class A domains protein 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use