


Product Name: Recombinant Human ERVK-10 Protein
Catalog #: PEH5834



Summary

Name	GAK10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P87889
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALEPFQTEED SISVSDAPGS CLIDCNENTR KKSQKETESL HCEYVAEPVM AQSTQNVDYN QLQEVYIPET LKLEGKGPEL MGPSESKPRG TSPLPAGQVL VRLQPQKQVK

image



Background

Alternative Names	Endogenous retrovirus group K member 10 Gag polyprotein (HERV-K10 Gag protein) (HERV-K107 Gag protein) (HERV-K_5q33.3 provirus ancestral Gag polyprotein) (Gag polyprotein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ERVK-10 Protein
Catalog #: PEH5834

