


Summary

Name	LEUTX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A8MZ59
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MHPSLATMGK LASKLQLDLS VVKIWFKNQR AKWKRRQQRQQ MQTRPSLGPA NQTTSVKKEE TPSAITTANI RPVSPGISDA

image



Background

Alternative Names	Paired-like homeodomain transcription factor LEUTX (Leucine-twenty homeobox) (Paired-like homeobox transcription factor LEUTX) (PRD-LIKE homeobox transcription factor LEUTX)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use