

Product Name: Recombinant Human KCNJ16 Protein
Catalog #: PEH5843



Summary

Name	KCJ16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPI9
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVEGTVRAQL LRYTEDSEGR MTMAFKDLKL VNDQIILVTP VTIVHEIDHE SPLYALDRKA VAKDNFEILV TFIYTG DSTG TSHQSRSSYV PREILWGHRF NDVLEVKKRY YKVNCLQFEG

image

SDS

Background

Alternative Names Inward rectifier potassium channel 16 (Inward rectifier K(+) channel Kir5.1)
(Potassium channel, inwardly rectifying subfamily J member 16)

Note

For Research Use Only , Not for Diagnostic Use