


Product Name: Recombinant Human ACSL5 Protein
Catalog #: PEH5854



Summary

Name	ACSL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9ULC5
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLDLNNQSVG IEGGARKGVS QKNNDLTSCC FSDAKTMYEV FQRGLAVSDN GPCLGYRKPN QPYRWLSYKQ VSDRAEYLGSL CLLHKGYKSS PDQFVGIFAQ NRPEWIISEL ACYTYSMVAV PLYDTLGPEA IVHIVNKADI AMVICDTPQK ALVLIGNVEK GFTPSLVKVII LMDPFDDDLK QRGEKSGIEI LSLYDAENLG KEHFRKPVPP SPEDLSVICF TSGTTGDPKG

image



Background

Alternative Names Long-chain-fatty-acid--CoA ligase 5 (EC 6.2.1.3) (Arachidonate--CoA ligase) (EC 6.2.1.15) (Long-chain acyl-CoA synthetase 5) (LACS 5)

Note

Product Name: Recombinant Human ACSL5 Protein
Catalog #: PEH5854



For Research Use Only , Not for Diagnostic Use