


Product Name: Recombinant Human ADAMTS3 Protein
Catalog #: PEH5860



Summary

Name	ATS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15072
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLIQEDEASP HILAIGNQAT GHYILNGKGE EAKSRTFIDL GVEWDYNIED DIESLHTDGP LHDPVIVLII

image



Background

Alternative Names	A disintegrin and metalloproteinase with thrombospondin motifs 3 (ADAM-TS 3) (ADAM-TS3) (ADAMTS-3) (EC 3.4.24.-) (Procollagen II N-proteinase) (PC II-NP) (Procollagen II amino propeptide-processing enzyme)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use