


Product Name: Recombinant Human DNAJC14 Protein
Catalog #: PEH5864



Summary

Name	DJC14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6Y2X3
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYSEGGNGSS SNFCHHCTSP ALGEDELEEE YDDEESLKFP SDFSRVSSGK KPPSRRQRHR FPTKEDTREG GRRDPRSPGR

image



Background

Alternative Names	DnaJ homolog subfamily C member 14 (DnaJ protein homolog 3) (Dopamine receptor-interacting protein of 78 kDa) (DRIP78) (Human DnaJ protein 3) (hDj-3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use