


Product Name: Recombinant Human FCGBP Protein
Catalog #: PEH5866



Summary

Name	FCGBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6R7
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQEVWADEL C QRRCTCNGAT HQVTCRDKQS CPAGERCSVQ NGLLG CYPDR FGTCQGSGDP HYVSFDGRRF DFMGTCTYLL VGSCGQNAAL PAFRVLVENE HRGSQTVSYT RAVRVEARGV KVAVRREYPG QVLVDDVLQY LPFQAADGQV QVFRQGRDAV VRTDFGLTVT YDWNARVTAK VPSSYAEALC GLCGNFNGDP ADDLALRGGG QAANALAFGN SWQEETRPGC GATEPGDCPK LDSLVAQQLQ SKNECGILAD PKGPFRECHS KLDPPQGA VRD CVYDRCLLPG QSGPLCDALA TYAAACQAAG ATVHPWRSEE LCPLSCPPHS HYEACSYGCP LSCGDLVPVG GCGSECHEGC VCDEGFALSG ESCLPLASCG CVHQGTYHPP GQTFYPGPGC

image



Background

Alternative Names	IgGFc-binding protein (Fcgamma-binding protein antigen) (FcgammaBP)
--------------------------	---

Product Name: Recombinant Human FCGBP Protein
Catalog #: PEH5866



Note

For Research Use Only , Not for Diagnostic Use