


Product Name: Recombinant Human FLT3 Protein
Catalog #: PEH5880



Summary

Name	FLT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P36888
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYFRKMENQD ALVCISEVP EPIVEWVLCD SQGESCKEES PAVVKKEEKV LHELFGTDIR CCARNELGRE CTRLFTIDLN QTPQTTLPL FLKVGEPPLWI RCKAVHVNHG FGLTWELENK ALEEGNYFEM STYSTNRTMI RILFAFVSSV ARNDTGYYTC SSSKHPSQSA LVTIVEKGFI NATNSSSEDYE

image



Background

Alternative Names	Receptor-type tyrosine-protein kinase FLT3 (EC 2.7.10.1) (FL cytokine receptor) (Fetal liver kinase-2) (FLK-2) (Fms-like tyrosine kinase 3) (FLT-3) (Stem cell tyrosine kinase 1) (STK-1) (CD antigen CD135)
--------------------------	--

Note

Product Name: Recombinant Human FLT3 Protein
Catalog #: PEH5880



For Research Use Only , Not for Diagnostic Use