


Product Name: Recombinant Human ARHGEF39 Protein
Catalog #: PEH5884



Summary

Name	ARG39_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N4T4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKRACTAREL LETERRYQEQ LGLVATYFLG ILKAKGTLRP PERQALFGSW ELIYGASQEL LPYLEGGCWG QGLEGFCRHL ELYNQFAANS ERSQTTLQEQ LKKNKGFRRF VRLQEGRPEF GGLQLQDLLP LPLQRLQQYE NLVVALAENT GPNSPDHQQQL TRAARLSET AQRVHTIGQK QKNDQHLRRV

image



Background

Alternative Names Rho guanine nucleotide exchange factor 39

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARHGEF39 Protein
Catalog #: PEH5884

