


Product Name: Recombinant Human CYP4F3 Protein
Catalog #: PEH5889



Summary

Name	CP4F3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08477
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IQERRRTLPS QGVDDFLQAK AKSKTLDFID VLLSKDEDEG KKLSDIEDIRA EADTFMFEGH DTTASGLSWV LYHLAKHPEY QERCRQEVQE

image



Background

Alternative Names	Cytochrome P450 4F3 (EC 1.14.14.1) (20-hydroxyeicosatetraenoic acid synthase) (20-HETE synthase) (CYPIVF3) (Cytochrome P450-LTB-omega) (Docosahexaenoic acid omega-hydroxylase CYP4F3) (EC 1.14.14.79) (Leukotriene-B(4) 20-monooxygenase 2) (Leukotriene-B(4) omega-hydroxylase 2) (EC 1.14.14.94)
--------------------------	---

Note

Product Name: Recombinant Human CYP4F3 Protein
Catalog #: PEH5889



For Research Use Only , Not for Diagnostic Use