


Product Name: Recombinant Human KLK7 Protein
Catalog #: PEH5893



Summary

Name	KLK7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P49862
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IIDGAPCARG SHPWQVALLS GNQLHCGGVL VNERWVLTAA HCKMNEYTVH LGSDTLGDRR AQRIKASKSF RHPGYSTQTH VNDLMLVKLN SQARLSSMVK KVRLPSRCEP PGTCTVSGW GTTSPDVTF PSDLMCVDVK LISPQDCTKV YKDLENSML CAGIPDSKKN ACNGDSGGPL VCRGTLQGLV SWGTFPCGQP NDPGVYTQVC KFTKWINDTM KKHR

image



Background

Alternative Names Kallikrein-7 (hK7) (EC 3.4.21.117) (Serine protease 6) (Stratum corneum chymotryptic enzyme) (hSCCE)

Note

Product Name: Recombinant Human KLK7 Protein
Catalog #: PEH5893



For Research Use Only , Not for Diagnostic Use