



Product Name: Recombinant Human KLK9 Protein
Catalog #: PEH5894



Summary

Name	KLK9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UKQ9
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGWADTRAIG AEECRPNSQP WQAGLFHLTR LFCGATLISD RWLLTAAHCR KPYLWVRLGE HHLWKWEGPE QLFRTDFFP HPGFNKDLSA NDHNDDIMLI RLPRQARLSP AVQPLNLSQT CVSPGMQCLI SGWGAVSSPK ALFPVTLQCA NISILENKLC HWAYPGHISD SMLCAGLWEG GRGSCQGDSG GPLVCNGTLA GVVSGGAEPC SRPRRPAVYT SVCHYLDWIIQ EIMEN

image



Background

Alternative Names Kallikrein-9 (EC 3.4.21.-) (Kallikrein-like protein 3) (KLK-L3)

Note

Product Name: Recombinant Human KLK9 Protein
Catalog #: PEH5894



For Research Use Only , Not for Diagnostic Use