


Product Name: Recombinant Human DAND5 Protein
Catalog #: PEH5909




Summary

Name	DAND5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N907
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLGMGRLQRG QDEVAAVTLP LNPQEVIQGM CKAVPFVQVF SRPGCSAIRL RNHLCFGHCS SLYIPGSDPT PLVLCNSCMP ARKRWAPVVL WCLTGSSASR RRVKISTMLI

image



SDS



WB

Background

Alternative Names	DAN domain family member 5 (Cerberus-like protein 2) (Cerl-2) (Cysteine knot superfamily 1, BMP antagonist 3) (Gremlin-3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use