


Product Name: Recombinant Human CA5A Protein
Catalog #: PEH5910



Summary

Name	CAH5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P35218
Host	E.coli
Species	Human
Molecular Mass	72 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WQTSNNTLHP LWTVPVSVPG GTRQSPINIQ WRDSVYDPQL KPLRVSYEAA SCLYIWNTGY LFQVEFDDAT EASGISGGPL ENHYRLKQFH FHWGAVNEGG SEHTVDGHAY PAELHLVHWN SVKYQNYKEA VVGENGLAVI GVFLKLGAAH QTLQRLVDIL PEIKHKDARA AMRPFDPSTL LPTCWDYWTY AGSLTTPPLT ESVTWIIQKE PVEVAPSQLS AFRTLFSAL GEEKMMVNN YRPLQPLMNR KVWASFQATN

image



Background

Alternative Names	Carbonic anhydrase 5A, mitochondrial (EC 4.2.1.1) (Carbonate dehydratase VA) (Carbonic anhydrase VA) (CA-VA)
--------------------------	--

Product Name: Recombinant Human CA5A Protein
Catalog #: PEH5910



Note

For Research Use Only , Not for Diagnostic Use