


Product Name: Recombinant Human CA5A Protein
Catalog #: PEH5911



Summary

| | |
|------------------------------|---|
| Name | CAH5A_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-SUMO & N-His |
| Accession # | P35218 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 42 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | WQTSNNTLHP LWTVPVSVPG GTRQSPINIQ WRDSVYDPQL KPLRVSYEAA SCLYIWNTGY LQVEFDDAT EASGISGGPL ENHYRLKQFH FHWGAVNEGG SEHTVDGHAY PAELHLVHWN SVKYQNYKEA VVGENGLAVI GVFLKLGAAH QTLQRLVDIL PEIKHKDARA AMRPFDPSTL LPTCWDYWTY AGSLTTPPLT ESVTWIIQKE PVEVAPSQLS AFRTLLFSAL GEEEKMMVNN YRPLQPLMNR KVVASFQATN |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Carbonic anhydrase 5A, mitochondrial (EC 4.2.1.1) (Carbonate dehydratase VA) (Carbonic anhydrase VA) (CA-VA) |
|--------------------------|--|

Product Name: Recombinant Human CA5A Protein
Catalog #: PEH5911



Note

For Research Use Only , Not for Diagnostic Use