


Product Name: Recombinant Human CDH12 Protein
Catalog #: PEH5926



Summary

Name	CAD12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55289
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FFVLEEYVGS EPQYVGKLHS DLDKGEGTVK YTLSGDGAGT VFTIDETTGD IHAIRSLDRE EKPFTYTLRAQ AVDIETRKPL EPESEFIKV QDINDNEPKF LDGPYVATVP EMSPVGAYVL QVKATDADDP TYGNSARVVY SILQGQPYFS IDPKTGIVRT ALPNMDREVK EQYQVLIQAK DMGGQLGGLA GTTIVNITLT

image



Background

Alternative Names	Cadherin-12 (Brain cadherin) (BR-cadherin) (Neural type cadherin 2) (N-cadherin 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CDH12 Protein
Catalog #: PEH5926

