


Product Name: Recombinant Human MMP11 Protein
Catalog #: PEH5932



Summary

Name	MMP11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P24347
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FVLSGGRWEK TDLTYRILRF PWQLVQEQVR QTMAEALKVW SDVTPLTFTE VHEGRADIMI DFARYWHGDD LPFDGPGGIL AHAFPPKTHR EGDVHFDYDE TWTIGDDQGT DLLQVAAHEF GHVLGLQHTT AAKALMSAFY TFRYPLSLSP DDCRGVQHLY GQPWPTVTSR TPALGPQAGI DTNEIAPLEP DAPPDACEAS FDAVSTIRGE LFFFKAGFVW RLRGGQLQPG YPALASRHWQ GLPSPVDAAF EDAQGHIWFF QGAQYWVYDG EKPVLGPAPL TELGLVRFPV HAALVWGPEK NKIYFFRGRD YWRFHPSTRR VDSPVPRRAT DWRGVPSEID AAFQDADGYA YFLRGRLYWK FDPVKVKALE GFPRLVGPDF FGCAEPANTF L

image



Background

Alternative Names Stromelysin-3 (SL-3) (ST3) (EC 3.4.24.-) (Matrix metalloproteinase-11)

Product Name: Recombinant Human MMP11 Protein
Catalog #: PEH5932



(MMP-11)

Note

For Research Use Only , Not for Diagnostic Use