


Product Name: Recombinant Human SLC16A4 Protein
Catalog #: PEH5933



Summary

Name	MOT5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15374
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSENNSGIKD KGSSLSAHGP EAHATETHCH ETEESTIKDS TTQKAGLPSK NLTVSQNNQSE EFYNGPNRNR LLLKSDEESD KVISWSCKQL

image



Background

Alternative Names	Monocarboxylate transporter 5 (MCT 5) (Monocarboxylate transporter 4) (MCT 4) (Solute carrier family 16 member 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use