


Product Name: Recombinant Human NUDT14 Protein
Catalog #: PEH5934



Summary

Name	NUDT14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST
Accession #	O95848
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MERIEGASVG RCAASPYLRP LTLHYRQNGA QKSWDFMKTH DSVTVLLFNS SRRSLVLVKQ FRPAVYAGEV ERRFPGSLAA VDQDGPRELQ PALPGSAGVT VELCAGLVDQ PGLSLEEVAC KEAWEECGYH LAPSDLRRVA TYWSGVGLTG SRQTMFYTEV TDAQRSGPGG GLVEEGELIE VVHLPLEGAQ AFADDPDIPK TLGVIFGVSW FLSQVAPNLD LQ

image



Background

Alternative Names	Uridine diphosphate glucose pyrophosphatase NUDT14 (UDPG pyrophosphatase) (UGPPase) (EC 3.6.1.45) (Nucleoside diphosphate-linked moiety X motif 14) (Nudix motif 14)
--------------------------	--

Product Name: Recombinant Human NUDT14 Protein
Catalog #: PEH5934



Note

For Research Use Only , Not for Diagnostic Use