


Product Name: Recombinant Human SPAM1 Protein
Catalog #: PEH5936



Summary

Name	HYALP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P38567
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGKFDEPLDM SLFSFIGSPR INATGQGVTI FYVDRLGYYP YIDSITGVTV NGGIPQKISL QDHLDAKAKD ITFYMPVDNL GMAVIDWEEW RPTWARNWKP KDVYKNRSIE LVQQQNVQLS LTEATEKAKQ EFEKAGKDFL VETIKLGKLL

image



Background

Alternative Names	Hyaluronidase PH-20 (Hyal-PH20) (EC 3.2.1.35) (Hyaluronoglucosaminidase PH-20) (Sperm adhesion molecule 1) (Sperm surface protein PH-20)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use