


Product Name: Recombinant Human MSL3P1 Protein
Catalog #: PEH5937



Summary

Name	MS3L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	POC860
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VARALSRRR YVCARDADAS RRRRRPFNYG LSIEEKEND ENSLSSSSDS SEDKDEKISE ECDIEEKTEV KEEPELQTKR

image



Background

Alternative Names	Putative male-specific lethal-3 protein-like 2 (MSL3-like 2) (Male-specific lethal-3 homolog 2) (Male-specific lethal-3 homolog pseudogene 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use