


Product Name: Recombinant Human PCGF2 Protein
Catalog #: PEH5946



Summary

Name	PCGF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35227
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YPLTEVPNGS NEDRGEVLEQ EKGALSDDEI VLSIEFYEG ARDRDEKKGP LENGDGDKEK TGVRFLRCPA AMTVMHLAKF LRNKMDVPSK YKVEVLYEDE PLKEYYTLMD IAYIYPWRRN GPLPLKYRVQ PACKRLTLAT VPTPSEGTNT SGASECESVS DKAPSPATLP

image



Background

Alternative Names Polycomb group RING finger protein 2 (DNA-binding protein Mel-18) (RING finger protein 110) (Zinc finger protein 144)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCGF2 Protein
Catalog #: PEH5946

