


Product Name: Recombinant Human CD1C Protein
Catalog #: PEH5949



Summary

Name	CD1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P29017
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASQEHVSFHV IQISFVNQS WARGQSGWL DELQTHGWDS ESGTIIFLHN WSKGNFSNEE LSDLELLFRF YLFGLTREIQ DHASQDYSKY PFEVQVKAGC ELHSGKSPEG

image



Background

Alternative Names T-cell surface glycoprotein CD1c (CD antigen CD1c)

Note

For Research Use Only , Not for Diagnostic Use