



Product Name: Recombinant Human MMP25 Protein
Catalog #: PEH5956



Summary

Name	MMP25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9NPA2
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MRPFYQGPVG DPDKYRLSQD DRDGLQQLYG KAPQTPYDKP TRKPLAPPPQ PPASPTHSPS FPIPDRCEGN FDAIANIRGE TFFFKGPWFW RLQPSGQLVS PRPARLHRFW EGLPAQVRVV QAAYARHRDG

image



Background

Alternative Names	Matrix metalloproteinase-25 (MMP-25) (EC 3.4.24.-) (Leukolysin) (Membrane-type matrix metalloproteinase 6) (MT-MMP 6) (MTMMP6) (Membrane-type-6 matrix metalloproteinase) (MT6-MMP) (MT6MMP)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MMP25 Protein
Catalog #: PEH5956

