


Product Name: Recombinant Human MMP28 Protein
Catalog #: PEH5958



Summary

Name	MMP28_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9H239
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQGRRPETQG PKYCHSSFDA ITVDRQQQLY IFKGSHFWEV AADGNVSEPR PLQERWVGLP PNIEAAAVSL NDGDFYFFKG GRCWRFRGPK PVWGLPQLCR AGGLPRHPDA ALFFPPLRRL ILFKGARYYV LARGGLQVEP YYPRSLQDWG GIPEEVSGAL PRPDGSIIF RDDRYWRLDQ AKLQATTSGR

image



Background

Alternative Names Matrix metalloproteinase-28 (MMP-28) (EC 3.4.24.-) (Epilysin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MMP28 Protein
Catalog #: PEH5958

