


Product Name: Recombinant Human KLK12 Protein
Catalog #: PEH5960



Summary

Name	KLK12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UKR0
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATPKIFNGTE CGRNSQPWQV GLFEGTSLRC GGVLDHRWV LTAAHCSGSR YWVRLGEHSL SQLDWTEQIR HSGFSVTHPG YLGASTSHEH DLRLRLRLP VRVTSSVQPL PLPNDCATAG TECHVSGWGI TNHPRNPFPD LLQCLNLSIV SHATCHGVYP GRITSNMVCA GGVPQGDAQ GDSGGPLVCG GVLQGLVSWG SVGPCGQDGI PGVYTYICKY VDWIRMIMRN N

image



Background

Alternative Names Kallikrein-12 (EC 3.4.21.-) (Kallikrein-like protein 5) (KLK-L5)

Note

Product Name: Recombinant Human KLK12 Protein
Catalog #: PEH5960



For Research Use Only , Not for Diagnostic Use