


Product Name: Recombinant Human KLK2 Protein
Catalog #: PEH5961



Summary

Name	KLK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20151
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CEKHSQPWQV AVYSHGWAHC GGVLVHPQWV LTAAHCLKKN SQVWLGRHNL FEPEDTGQRV PVSHSFPHPL YNMSLLKHQS LRPDEDSSHD LMLLRLSEPA KITDVVKVLG LPTQEPALGT TCYASGWGSI EPEEFLRPRS LQCVSLHLLS NDMCARAYSE KVTEFMLCAG LWTGGKDTG GDSGGPLVCN GVLQGITSWG PEPCALPEKP AVYTKVVHYR

image



Background

Alternative Names Kallikrein-2 (EC 3.4.21.35) (Glandular kallikrein-1) (hGK-1) (Tissue kallikrein-2)

Note

Product Name: Recombinant Human KLK2 Protein
Catalog #: PEH5961



For Research Use Only , Not for Diagnostic Use