


Product Name: Recombinant Human KLK4 Protein
Catalog #: PEH5969



Summary

Name	KLK4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5K2
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IINGEDCSPH SQPWQAALVM ENELFCSGVL VHPQWVLSAA HCFQNSYTIG LGLHSLEADQ EPGSQMVEAS LSVRHPEYNR PLLANDLMLI KLDESVSSED TIRSISIASQ CPTAGNSCLV SGWGLLANGR MPTVLQCVNV SVVSEEVCSK LYDPLYHPSM FCAGGGHDQK DSCNGDSGGP LICNGYLQGL VSFGKAPCGQ VGVPGVYTNL CKFTEWIEKT VQAS

image



Background

Alternative Names	Kallikrein-4 (EC 3.4.21.-) (Enamel matrix serine proteinase 1) (Kallikrein-like protein 1) (KLK-L1) (Protease) (Serine protease 17)
--------------------------	---

Note

Product Name: Recombinant Human KLK4 Protein
Catalog #: PEH5969



For Research Use Only , Not for Diagnostic Use