


Product Name: Recombinant Human PARP16 Protein
Catalog #: PEH5977



Summary

Name	PAR16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5Y8
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAAGRDMMLAA DLRCSLFASA LQSYKRDSVL RPFASYARG DCKDFEALLA DASKLPNLKE LLQSSGDNHK RAWDLVSWIL SSKVLTIHSA GKAEFEKIQK LTGAPHTPVP APDFLFEIEY FDPANAKFYE TKGERDLIYA FHGSRLNFH SIIHNGHLCH LNKTSLFEGEG TYLTSDLALA LIYSPHGHGW QHSLGPILS CVAVCEVIDH PDVKCQTKKK DSKEIDRRRA RIKHSEGGDI PPKYFVVTNN QLLRVKYLLV YSQKPPKRAS

image



Background

Alternative Names	Protein mono-ADP-ribosyltransferase PARP16 (EC 2.4.2.-) (ADP-ribosyltransferase diphtheria toxin-like 15) (Poly [ADP-ribose] polymerase 16) (PARP-16)
--------------------------	---

Product Name: Recombinant Human PARP16 Protein
Catalog #: PEH5977



Note

For Research Use Only , Not for Diagnostic Use