


Product Name: Recombinant Human EID2B Protein
Catalog #: PEH5983



Summary

Name	EID2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96D98
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SELPGDSSVP QVG TASGVSD VLRGAVGGGV RVQEAREGPV AEAARSMARM PGPVPGPIPS SVPGLASAPD PHQQLAFLEI NRQLLFREYL DGSSMIPVRL LRDFEERRRL FVEGCKAREA AFDADPPQMD FAAVAFTVAL

image



Background

Alternative Names	EP300-interacting inhibitor of differentiation 2B (EID-2B) (EID-2-like inhibitor of differentiation 3) (EID-3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use