


Product Name: Recombinant Human GPX5 Protein
Catalog #: PEH5985



Summary

Name	GPX5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75715
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPKQEKMKMD CHKDEKGTIY DYEAIALNKN EYVSFKQYVG KHILFVNVAT YCGLTAQYPE LNALQEELKP YGLVVLGFPC NQFGKQEPGD NKEILPGLKY VRPGGGFVPS FQLFEKGDVN GEKEQKVFSF LKHSCPHPSE ILGTFKSISW DPVKVHDIRW NFEKFLVGPD GIPVMRWSHR ATVSSVKTDI LAYLKQFKTK

image



Background

Alternative Names	Epididymal secretory glutathione peroxidase (EC 1.11.1.9) (Epididymis-specific glutathione peroxidase-like protein) (EGLP) (Glutathione peroxidase 5) (GPx-5) (GSHPx-5)
--------------------------	---

Note

Product Name: Recombinant Human GPX5 Protein
Catalog #: PEH5985



For Research Use Only , Not for Diagnostic Use