


Product Name: Recombinant Human GPX6 Protein
Catalog #: PEH5986



Summary

Name	GPX6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P59796
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQTLPQNRK VDCNKGVTGT IYEGALTIN GEEYIQFKQF AGKHVLFVNV AAYUGLAAQY PELNALQEEL KNFGVIVLAF PCNQFGKQEP GTNSEILLGL KYVCPGSGFV PSFQLFEKGD VNGEKEQKVF TFLKNSCPPT SDLLGSSSQL FWEPMKVHDI RWNFEKFLVG PDGVPVMHWF HQAPVSTVKS DILEYLKQFN TH

image



Background

Alternative Names Glutathione peroxidase 6 (GPx-6) (GSHPx-6) (EC 1.11.1.9)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GPX6 Protein
Catalog #: PEH5986

