


Product Name: Recombinant Human GPX7 Protein
Catalog #: PEH5987



Summary

Name	GPX7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SL4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQEQDFYDFK AVNIRGKLVS LEKYRGSVSL VVNVASECGF TDQHYRALQQ LQRDLGPHHF NVLAFPCNQF GQQEPDSNKE IESFARRTYS VSFPMFSKIA VTGTGAHPAF KYLAQTSGKE PTWNFWKYL V APDGKVVGAW DPTVSVEEVR PQITALVRKL ILLKREDL

image



Background

Alternative Names Glutathione peroxidase 7 (GPx-7) (GSHPx-7) (EC 1.11.1.9) (CL683)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GPX7 Protein
Catalog #: PEH5987

